



Membrane filters millipore S-pak type HA, 47 mm, 0.45 um pore size, pk/600

Številka izdelka: 2936500

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

For membrane filtration in microbiological analysis of aqueous samples. "HA" signifies certified compliance with Standard Methods criteria with pore size of 0.45 um. White, cellulose and mixed-ester composition. S-Pak (filter only) is sterile. Diameter: 4

Specifikacije

Diameter:	47 mm mm
Material:	Mixed-Cellulose Ester
Prostornina / velikost paketa:	600
Velikost por:	0.45 µm