



Fixed volume, Micropipette, 1000 µL

Številka izdelka: 2946010

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and color-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Specifikacije

Kapaciteta:	1000 µL
Natančnost:	0,5 %
Prostornina / velikost paketa:	1 each