



QL filters, pore size 8 - 12 μm , diameter 185 mm, circular, 100/pk, qualitative

Številka izdelka: 50659
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Specifikacije

Diameter:	185 mm
Material:	Cellulose
Prostornina / velikost paketa:	100 /pk