



Cylinder, graduated, 25 mL +- 0.3 mL, 0.5 mL divisions (white markings)

Številka izdelka: 50840
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Specifikacije

Certificiranje:	Class B
Kapaciteta:	25 mL
Material ohišja:	Glass
Merilna metoda:	To Contain
Prostornina / velikost paketa:	1 each
Vrsta:	eno metričen