



## **Buchner Funnel 320 mL approx. volume, Porcelain**

|                              |         |
|------------------------------|---------|
| Številka izdelka:            | 55087   |
| EUR Cena:                    | Kontakt |
| Datum pošiljanja ni prikazan |         |

Durable porcelain construction is tough, inert and highly chemical resistant. Molded-in perforated filter plate. Fully autoclavable.

---

### **Specifikacije**

|                     |           |
|---------------------|-----------|
| Approximate Volume: | 320 mL    |
| Kapaciteta:         | 320 mL    |
| Material ohišja:    | Porcelain |