



Titration Application Pack Water: pH, Alkalinity & Hardness (ISE)

Številka izdelka: **AP0005.AT1222**

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Application Pack for the determination of pH, Alkalinity and Calcium & Magnesium Hardness in water. Suited for Titrator Model AT1000 Series, Automatic Titrator Model AT1222.

Specifikacije

Aplikacija:	Municipal Applications
Garancija:	24 mesecev
Parameter:	pH/Alkalinity & Hardness
Vrsta:	Application Pack
Vsebina paketa:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.
Združljivi instrumenti:	AT1222

Vsebina paketa

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.