



Titration Application Pack Biogas: FOS/TAK

Številka izdelka: **AP0006.AT1102**

EUR Cena: Kontakt

Odprema v 2 tednih

Application Pack for the determination of FOS/TAK (Biogas). Suited for Titrallab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

Specifikacije

Aplikacija:	Environmental Applications
Garancija:	24 mesecev
Parameter:	FOS/TAC (Biogas)
Vrsta:	Application Pack
Vsebina paketa:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.
Združljivi instrumenti:	AT1102
	AT1112
	AT1122
	AT1222

Vsebina paketa

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.