



Titration Application Pack Food & Beverage: Salt

Številka izdelka: **AP0010.AT1112**

EUR Cena: Kontakt

Odprema v 2 tednih

Application Pack for the determination of Salt in food products. Suited for Titrallab AT1000 Series, Automatic Titrator Models AT1112/AT1122/AT1222.

Specifikacije

Aplikacija:	Food & Beverage Applications
Garancija:	24 mesecev
Parameter:	Salt in food products
Vrsta:	Application Pack
Vsebina paketa:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.
Združljivi instrumenti:	AT1112
	AT1122
	AT1222

Vsebina paketa

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.