



Titration Application Pack Petrochemistry: TAN (Total Acid Number)

Številka izdelka: **AP0015.AT1102**

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Application Pack for the determination of the Total Acid Number in Petrochem. Suited for Titrilab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

Specifikacije

Aplikacija:	Petrochemical Applications
Garancija:	24 mesecev
Parameter:	TAN (Total Acid Number)
Vrsta:	Application Pack
Vsebina paketa:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.
Združljivi instrumenti:	AT1102
	AT1112
	AT1122
	AT1222

Vsebina paketa

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.