



Buffer Solution Kit, Colour-coded, pH 4.01, pH 7.00 and pH 10.01, 500 mL

Številka izdelka: 2947600
EUR Cena: Kontakt
Na voljo

Save when you purchase pH buffer 3-packs!

pH Buffer Standard Solution for pH electrode calibration. Pack of 3 Colour-coded pH buffers. pH is ± 0.02 pH units. 500 mL bottles of pH 4.01, 7.00 & 10.01

33% Cost Reduction over single bottle purchase

Specifikacije

Barva kodirana?:	Da
Koncentracija:	4,01 pH - 10,01 pH
Opis:	pH Buffer Calibration reagent
Package settings:	Bottle
pH vrednost(i):	4.01, 7.00 & 10.01 ± 0.02 pH at 25°C
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Prostornina / velikost paketa:	3 pk of 500 mL
Sledljivo po NIST-u:	Da