



PhosVer 3 Phosphate Reagent Powder Pillows 0.02-2.50mg/L PO₄

Številka izdelka:

212599

Nearno

EUR Cena:

Kontakt

Na voljo



Powder reagent in sealed package for the determination of Phosphorus, reactive (25mL sample; 100 tests).

Specifikacije

Ime metode:	Ascorbic Acid
Inštrument:	DR2000, DR3000, DR900
Merilno območje:	0.02 - 2.50 PO ₄
Parameter:	Phosphorus
Platforma:	Praškasti reagenti
Prostornina / velikost paketa:	100 /pk
Številka metode:	8048
Število testov:	100