



Nitrate standard, 100 mg/L NO₃ (22.6 mg/L NO₃-N), 500 mL

Številka izdelka:	LCW826
EUR Cena:	Kontakt
Na voljo	

Standard solution for online nitrate determination

Specifikacije

Inštrument:	Nitratax plus sc, Nitratax clear sc, NT3100sc, NT3200sc
Koncentracija:	100 mg/L NO ₃ (22.6 mg/L NO ₃ -N)
Parameter:	Nitrate
Platforma:	Raztopina
Prostornina / velikost paketa:	500 mL