



Nitrate standard, 400 mg/L NO₃ (90.4 mg/L NO₃-N), 500 mL

Številka izdelka: LCW863

EUR Cena: Kontakt

Na voljo

Specifikacije

Inštrument:	Nitratax plus sc, Nitratax clear sc, NT3100sc, NT3200sc
Koncentracija:	400 mg/L NO ₃ (90.4 mg/L NO ₃ -N)
Parameter:	Nitrate
Platforma:	Raztopina
Prostornina / velikost paketa:	500 mL