



Pole mounting hardware Filterprobe/Filtrax, 24 cm bracket, SS pole 2 m

Številka izdelka:

LZY714.99.42020

EUR Cena:

Kontakt

Na voljo

Specifikacije

Diameter:	40 mm
Dolžina:	Pole: 2000 mm
Inštrument:	suitable for: FILTERPROBE sc; FILTRAX; SIGMATAX . non verbal instruction sheet
Jeziki v navodilih za uporabo:	
Material:	Wetted parts: pole SS 1.4301; seals EPDM; screw set SS.
Teža:	Lifting part only: 5.5 kg