



Iron Reagent Powder Pills, 0.02-3.00mg/L Fe, FerroVer

Številka izdelka:
EUR Cena:
Odprema v 2 tednih

85428
Kontakt

Nearno



Powder reagent in sealed package, for 25mL sample.

FerroVer Iron Reagent: For the determination of total soluble Iron including all ferric (Fe_3^+) and ferrous (Fe_2^+) Iron as well as complexed Iron compounds (Fe-EDTA). Perfect for many applications and a wide range of concentrations. Determines Iron in the presence of Copper.

Specifikacije

EPA compliant:	Da
Ime metode:	FerroVer
Merilno območje:	0,02 - 3,00 mg/L Fe
Platforma:	Praškasti reagenti
Prostornina / velikost paketa:	1000 /pk
Številka metode:	8008
Število testov:	1000