



Reagent powder pillows, FerroVer iron, 10 mL, pk/1000

Številka izdelka:

2105728

EUR Cena:

Kontakt

Na voljo

Nearno



Powder reagent in sealed package, pk/1000. For 10mL sample.

Specifikacije

EPA compliant:	Da
Ime metode:	FerroVer
Inštrument:	All except DR/3000, DR/2000
Merilno območje:	0,02 - 3,00 mg/L Fe
Parameter:	Iron
Platforma:	Prašasti reagenti
Pogoji skladiščenja:	10 °C - 25 °C
Prostornina / velikost paketa:	1000 /pk
Številka metode:	8008
Število testov:	1000