



MolyVer 1 Molybdenum Reagent Powder Pillows, pk/100

Številka izdelka:

2604299

Nearno

EUR Cena:

Kontakt

Odprema v 2 tednih



Powder reagent in sealed packages for the determination of Molybdenum, 0.3-40.0mg/L Mo. For 10mL sample.

Specifikacije

Ime metode:	Mercaptoacetic Acid
Parameter:	Molybdenum
Platforma:	Prašasti reagenti
Prostornina / velikost paketa:	100 /pk
Številka metode:	8036
Število testov:	100