



Sodium periodate powder pillows for manganese, 25 mL, pk/100

Številka izdelka:

98499

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Powder reagent in sealed package for the determination of Manganese, cold periodate method, and Iron titrations. For 25 ml sample.

For determination of manganese by the periodate oxidation method. EPA approved for reporting wastewater analyses (digestion required). Hach method 8034. Sample size: 25 mL Pack of 100 Powder Pillows.

Specifikacije

Ime metode:	Periodate Oxidation
Platforma:	Praškasti reagenti
Prostornina / velikost paketa:	100 /pk
Število testov:	100