



QAC reagent 2, 0.2- 5.0 mg/L CTAB, 25 powder pillows

Številka izdelka:

2401268

EUR Cena:

Kontakt

Odprema v 2 tednih

Nearno



Powder reagent in sealed package for the determination of Quaternary Ammonium Compounds (QAC).

Specifikacije

Ime metode:	Direct Binary Complex
Parameter:	Quaternary Ammonium Compounds
Platforma:	Prašasti reagenti
Prostornina / velikost paketa:	25/pk
Številka metode:	8337
Število testov:	12.5