



Silver nitrate digital titrator cartridge, 0.2256 N, 1 piece

Številka izdelka:

1439601

EUR Cena:

Kontakt

Na voljo

Nearno



Ready-to-use titrant solution for the analysis of Chloride (Silver Nitrate Method), 10 - 100 mg/L Cl⁻.

Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Ime metode:	Silver Nitrate
Merilno območje:	10 - 100 mg/L Cl ⁻
Parameter:	Chloride
Platforma:	Digitalni Titrator
Številka metode:	8207
Število testov:	40 - 100
Titrant:	0.2256 N AgNO ₃