



## Silver nitrate digital titrator cartridge, 1.128 N, 1 piece

Številka izdelka:

1439701

EUR Cena:

Kontakt

Na voljo

**Nearno**



Ready-to-use titrant solution for the analysis of Chloride (Silver Nitrate Method), 10 - 10000 mg/L Cl<sup>-</sup>.

For chloride determination by Digital Titrator argentometric titration. Hach method 8207. Digits = mg/L Cl<sup>-</sup> in 50 mL sample. Approximately 40-100 tests

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Specifikacije

|                  |                                  |
|------------------|----------------------------------|
| Ime metode:      | Silver Nitrate                   |
| Merilno območje: | 100 - 10000 mg/L Cl <sup>-</sup> |
| Parameter:       | Chloride                         |
| Platforma:       | Digitalni Titrator               |
| Številka metode: | 8207                             |
| Število testov:  | 40 - 100                         |
| Titrant:         | 1.128 N AgNO <sub>3</sub>        |