



Sodium Hydroxide Digital Titrator Cartridge, 0.3636 N

Številka izdelka: 1437801

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Nearno



Ready-to-use titrant solution for the determination of Carbon Dioxide, 10-100mg/l CO₂.

For carbon dioxide determination by Digital Titration

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Ime metode:	Phenolphthalein
Merilno območje:	10 - 100 mg/L CO ₂
Parameter:	Carbon Dioxide
Platforma:	Digitalni Titrator
Številka metode:	8205
Število testov:	40 - 100
Titrant:	0.3636 N NaOH