



Digital titrator cartridge, sulfuric acid, 1.600 N, 1 piece

Številka izdelka:

1438901

EUR Cena:

Kontakt

Na voljo

Nearno



Ready-to-use titrant solution for acid-base and alkalinity titrations.

Specifikacije

EPA compliant:	Da
Ime metode:	Verde de bromcresol/rojo de metilo
Merilno območje:	100 - 4000 mg/L CaCO ₃
Method Number 1 Chemistry:	8203
Parameter:	Alkalinity, Total
Platforma:	Digitalni Titrator
Število testov:	40 - 100
Titant:	1.600 N H ₂ SO ₄