



Alkalinity Reagent Set for Digital Titrator

Številka izdelka:

2271900

EUR Cena:

Kontakt

Odprema v 2 tednih

Nearno



Ready to use titrants in easily exchangeable, resealable cartridges. For reliable results with the Digital Titrator in the field and in the laboratory.

For alkalinity determination by sulfuric acid titration with the Digital Titrator. Hach method 8203. Range: 10-4,000 mg/L CaCO₃.

Set contains Phenolphthalein Indicator Powder Pillows, Bromcresol Green-Methyl Red Indicator Powder Pillows and Sulfuric Acid Digital Titrator Cartridges, 0.1600 N and 1.600 N. Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Ime metode:	Phenolphthalein and Total Alkalinity
Merilno območje:	10 - 4000 mg/L CaCO ₃
Opis:	Alkalinity Reagent Set for Digital Titrator
Parameter:	Alkalinity
Platforma:	Digitalni Titrator
Posebni zapiski:	Set contains: Bromcresol Green Methyl Red Indicator Powder Pillows (94399), Phenolphthalein Indicator (94299), Sulfuric Acid Titration Cartridge 1.600 N (1438901), Sulfuric Acid Titration Cartridge 0.1600 N (1438801)
Številka metode:	8203
Število testov:	100