



Hexamethylenetetramine Reagent Powder Pillows, pk/50

Številka izdelka:

187866

Nearno

EUR Cena:

Kontakt

Na voljo



Powder reagent in sealed package for the determination of Aluminium, ECR Method, 0-0.220 mg/L Al.

Specifikacije

Ime metode:	Eriochrome Cyanine R
Inštrument:	DR2700/DR2800/DR3800/DR3900/DR5000/DR6000
Parameter:	Aluminum
Platforma:	Prašasti reagenti
Prostornina / velikost paketa:	50 /pk
Številka metode:	8326