



Porphyrin 1 Reagent Powder Pillows

Številka izdelka:	2603549
EUR Cena:	Kontakt
Na voljo	

Powder reagent in sealed package for the determination of Copper with Porphyrin method, 2-210µg/L Cu.

Specifikacije

Ime metode:	Porphyrin
Parameter:	Copper
Platforma:	Praškasti reagenti
Prostornina / velikost paketa:	100 /pk
Številka metode:	8143
Število testov:	50 - 100

Obvezni pripomočki

- Copper Masking Reagent Powder Pillows, 10 mL, pk/100 (Item 2603449)
- Porphyrin 2 Reagent Powder Pillows (Item 2603649)