



Porphyrin 2 Reagent Powder Pillows

Številka izdelka:

2603649

Nearno

EUR Cena:

Kontakt

Na voljo



Powder reagent in sealed package for the determination of Copper with Porphyrin method, 2-210 µg/l Cu.

For Copper determination by the Porphyrin method. Hach Method 8143 Range: 2 - 210 µg/L Cu. Sample Size: 10 mL Pack of 100 Powder Pillows.

Specifikacije

Ime metode:	Porphyrin
Platforma:	Praškasti reagenti
Prostornina / velikost paketa:	100 /pk
Število testov:	50 - 100