



Water, deionised and demineralised, 4L

Številka izdelka:	27256
EUR Cena:	Kontakt
Na voljo	

Demineralised Type 2 Water

Specific Conductance: $\leq 1 \mu\text{S/cm}$,

Silica (SiO_2): $\leq 30 \mu\text{g/L}$,

Chloride (Cl^-): $\leq 20 \mu\text{g/L}$,

Hardness: Negative,

Consumption of KMnO_4 : Negative

Specifikacije

Platforma:	Raztopina
Prostornina / velikost paketa:	4 L