



Be Right™



Praškesti reagenti za prosti klor DPD, 10 mL, paket/1000

Številka izdelka: 2105528

EUR Cena: Kontakt

Odprema v 1 tednu

Nearno



Praškesti reagenti za prosti klor DPD v blazinicah za določanje prostega klora. Razpon 0,02–2,00 mg/L Cl₂. Za 10-mL vzorec.

Specifikacije

EPA compliant:	Da
Ime metode:	DPD
Inštrument:	Vsi razen DR/3000, DR/2000
Komentar:	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Merilno območje:	0,02 - 2,00 mg/L Cl ₂
Parameter:	Klor, prosti, klorov dioksid.
Platforma:	Praškesti reagenti
Prostornina / velikost paketa:	1000 /pk
Številka metode:	8021
Število testov:	1000
Volume / package size e Post Text:	/pk