



Sodium Thiosulfate Digital Titrator Cartridge, 0.113 N

Številka izdelka: 2267301
EUR Cena: Kontakt
Na voljo

Ready-to-use titrant solution for the analysis of Total Chlorine, 20-2000 mg/l Cl₂.

For chlorine determination by iodometric Digital Titrator titration. 0.113 N, Digits x 0.1 = mg/L Cl₂ in a 5 mL sample, Digits x 0.000141 = meq. Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Ime metode:	Iodometric
Merilno območje:	20 - 2000 mg/L Cl ₂
Parameter:	Chlorine, Total
Platforma:	Digitalni Titrator
Številka metode:	8209
Število testov:	40 - 100
Titrant:	0.113 N Na ₂ S ₂ O ₃