



Chlorine (Free and Total) Reagent Set, 0-3.5 mg/L Cl₂

Številka izdelka:

2438800

Nearno

EUR Cena:

Kontakt

Odprema v 2 tednih



For determination of Free and Total Chlorine by the DPD method, using visual test kits. Set contains: DPD Free Chlorine PP, pk/100 and DPD Total Chlorine PP pk/100. Range: 0 - 3.5 mg/L Cl₂. Sample Size: 5 mL.

Boric Acid free!

Specifikacije

Ime metode:	DPD
Merilno območje:	5 mL
Število testov:	100