



DPD Free Chlorine Reagent Powder Pillows, 5 mL, pk/1000

Številka izdelka:

1407728

Nearno

EUR Cena:

Kontakt

Na voljo



For determination of free chlorine by the DPD method.

For use with Hach's Chlorine Visual Test Kits.

Sample Size: 5 mL. Range: up to 3.5 mg/L Cl₂.

Pack of 1000 powder pillows.

Boric Acid free!

Specifikacije

Ime metode:	DPD
Merilno območje:	5 mL
Platforma:	Prašasti reagenti
Pogoji skladiščenja:	10 - 25 °C
Število testov:	1000