



DPD Free Chlorine Reagent Powder Pillows, 5 mL, pk/50

Številka izdelka:

2197846

Nearno

EUR Cena:

Kontakt

Na voljo



For determination of free chlorine with Hach's Chlorine Visual Test Kits- 2060300. Sample Size: 5 mL. Range: up to 2.5 mg/L Cl_2 . Pack of 50 Powder Pillows.

Boric Acid free!

Specifikacije

Ime metode:	DPD
Merilno območje:	5 mL
Parameter:	Chlorine, Free
Platforma:	Praškasti reagenti
Pogoji skladiščenja:	10 - 25 °C
Število testov:	50