



Chlorophosphonazo Indicator Solution Pillows, pk/100

Številka izdelka:

2589599

Nearno

EUR Cena:

Kontakt

Odprema v 2 tednih



For the determination of Ultra Low Range Hardness, 8-1000 µg/L CaCO₃. Chlorophosphonazo indicator solution is packaged in unit-dose solution pillows.

For ultra low range hardness determination by the Chlorophosphonazo method. Hach Method 8374. Range: 4 - 1000 µg/L Ca & Mg as CaCO₃. Chlorophosphonazo solution reagent packaged in solution pillows. Requires CDTA Reagent, Cat. No. 2589636. 100 Tests.

Specifikacije

Ime metode:	Chlorophosphonazo
Method Number 1 Chemistry:	8374
Parameter:	Calcium and Magnesium Hardness
Prostornina / velikost paketa:	100 /pk
Število testov:	100