



Be Right™



Reagent set, hardness (total), HR, digital titrator

Številka izdelka: 2448100

EUR Cena: Kontakt

Na voljo

Nearno



For the determination of Total Hardness, 100-4000 mg/l CaCO₃, EDTA Method.

Reagent set for the determination of Total Hardness by EDTA Digital Titrator titration (Hach method 8213). Range: 100-4,000 mg/L as CaCO₃. Set contains: Buffer Solution (Hardness 1), ManVer® 2 Hardness Indicator Powder Pillows and 0.800 M EDTA Digital Titrator Cartridge for 40 - 100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Ime metode:	EDTA Titration
Merilno območje:	100 - 4000 mg/L CaCO ₃
Platforma:	Digitalni Titrator
Številka metode:	8213
Število testov:	100