



Phenol 2 Reagent Powder Pillows

Številka izdelka:	183699
EUR Cena:	Kontakt
Na voljo	

Powder reagent in sealed packages for the determination of phenols by the 4-Aminoantipyrine method. pk/100

Specifikacije

Ime metode:	4-Aminoantipyrine
Opis:	Phenol 2 Reagent Powder Pillows
Parameter:	Phenols
Platforma:	Praškasti reagenti
Pogoji skladiščenja:	10 °C - 25 °C
Prostornina / velikost paketa:	100
Številka metode:	8047
Število testov:	50