



Chloride Reagent Set, 10-10000mg/L Cl

Številka izdelka:

2288000

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



For Silver Nitrate Digital Titrator Titration.

For determination of chloride by Silver Nitrate Digital Titrator Titration. Hach Method 8207. Ranges: 10-100 and 100-10,000 mg/L as Cl⁻. Set contains: Chloride 2 Indicator Powder Pillows and two Silver Nitrate Titration Cartridges (0.2256 N and 1.128 N). Approximately 40 to 100 tests per range.

Specifikacije

Ime metode:	Silver Nitrate Titration
Merilno območje:	10 - 10000 mg/L Cl ⁻
Platforma:	Digitalni Titrator
Številka metode:	8207
Število testov:	100