



Standardna raztopina amonijaka, 100 mg/L NH₃-N, 500-mililitrska

Številka izdelka: 2406549

EUR Cena: Kontakt

Na voljo

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

Enoparametrski standard za zagotavljanje analizne kakovosti.

Za preverjanje natančnosti pri določanju amonijevega dušika.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Koncentracija:	100 mg/L ±2 mg/L NH ₃ -N
Opis:	Dušik – amonijak, standardna raztopina, 100 mg/L NH ₃ -N, 500-mililitrska
Package settings:	Bottle
Parameter:	Amonijev dušik
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	500 mL 500