



Stablcal Turbidity Standards Calibration Kit for 2100P

Številka izdelka:

2659405

Nearno

EUR Cena:

Kontakt

Na voljo



Four ready to use, stabilised Formazin Primary standard solutions. Kit includes sealed vials of each: <0.1 NTU, 20 NTU, 100 NTU, 800 NTU .

Specifikacije

Inštrument:	2100P
Koncentracija:	< 0,1, 20, 100, 800 NTU
Parameter:	Turbidity
Prostornina / velikost paketa:	Set of 4 sealed vials