



Calibration set, Stablcal turbidity standards, 2100P portable turbidimeter, 100 mL bottles

Številka izdelka:

2659410

Nearno

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan



Four ready to use, stabilised Formazine Primary standard solutions in bottles. Kit includes 100 mL of each: <0.1 NTU, 20 NTU, 100 NTU, 800 NTU .

Specifikacije

Koncentracija:

< 0,1, 20, 100, 800 NTU

Prostornina / velikost paketa:

set of 4