



Stablcal Turbidity Standards Calibration Kit, 0 - 7500 NTU

Številka izdelka:
EUR Cena:
Odprema v 2 tednih

2659505
Kontakt

Nearno



Six ready to use, stabilised Formazine Primary standard solutions. Kit includes sealed vials of each: <0.1 NTU, 20 NTU, 200 NTU, 1000 NTU, 4000 NTU, 7500 NTU. Calibration Kit, for 2100AN, 2100AN IS, TL2350 and TL2360 IS Laboratory Turbidimeter.

Specifikacije

Koncentracija:	< 0,1, 20, 200, 1000, 4000, 7500 NTU
Prostornina / velikost paketa:	6 /pk