



Calibration kit, Stablcal turbidity standards, 2100AN / AN IS turbidimeter, 100 mL bottles

Številka izdelka:

2659510

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Ready for use, stabilised Formazine Primary Standards in bottles, Set: <0.1, 20, 200, 1000, 4000, 7500 NTU

Specifikacije

Inštrument:	2100 AN / AN IS
Koncentracija:	< 0,1, 20, 200, 1000, 4000, 7500 NTU
Opis:	Stablcal Calibration Kit for 2100 AN / AN IS (100 mL Bottles)
Parameter:	Turbidity
Platforma:	Raztopina
Pogoji skladiščenja:	5 °C - 25 °C
Prostornina / velikost paketa:	Set of 6