



Turbidity Standard, Stablcal, 800 NTU, 500 mL

Številka izdelka:

2660549

Nearno

EUR Cena:

Kontakt

Odprema v 2 tednih



Ready to use, stabilised Formazine Primary standard solution in 500 mL bottle, 800 NTU.

StablCal® Turbidity Standard for all Turbidimeters. USEPA accepted for calibration purposes. StablCal® Stabilized Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. NTU Value: 800. Bottle: 500 mL.

Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares StablCal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

Based on Formazin Primary Standards

StablCal Stabilized Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. StablCal is a stabilized version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. StablCal standards are based on decades of experience in turbidity measurement.

Specifikacije

Koncentracija:	800 NTU
Platforma:	Raztopina
Prostornina / velikost paketa:	500 mL