



Secondary Turbidity Standards Kit for 2100 Portable Turbidimeters

Številka izdelka:

2464105

Nearno

EUR Cena:

Kontakt

Na voljo



Kit includes 1 sealed vial each, 0-10 NTU, 10-100 NTU, 100-1000 NTU.

Gelex Secondary Standards Kit, for 2100P and 2100Q Portable Turbidimeters (includes 1 sealed vial each, 0-10 NTU, 10-100 NTU, 100-1000 NTU).

Specifikacije

Prostornina / velikost paketa:

Set of 3