



Pufrska raztopina za pH; pH 7,00 (500 mL)

Številka izdelka: 1222249
EUR Cena: Kontakt
Na voljo

pH Buffer standard solutions for laboratory and field analyses

Kakovostna pufrska raztopina; brezbarvna. Sledljivost po NIST. pH pri 25 °C je $7,00 \pm 0,02$.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Barva kodirana?:	Ne
Koncentracija:	7,00 pH
Obseg:	500 mL
Opis:	Pufrska raztopina (brezbarvna) (NIST)
Package settings:	Bottle
pH vrednost(i):	7.00
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf-life after opening, following recommended storage and use: 6 months
Prostornina / velikost paketa:	500 mL
Sledljivo po NIST-u:	Da