



Standardna raztopina natrijevega klorida, 1990 µS/cm

Številka izdelka: 210542
EUR Cena: Kontakt
Na voljo

Conductivity standard solutions for laboratory and field analyses

995 mg/L TDS, 1000 mg/L NaCl. Za sonde za merjenje prevodnosti.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Formula:	NaCl
Koncentracija:	NA
Opis:	Conductivity calibration solution
Package settings:	Bottle
Parameter:	Conductivity
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	100 mL
Sledljivo po NIST-u:	Da