



Buffer Powder Pillows, pH 7.00 Colour-coded Yellow, pk/250

Številka izdelka:

2227064

EUR Cena:

Kontakt

Odprema v 2 tednih

Nearno



pH Buffer standard solutions for laboratory and field analyses

pH Buffer Powder Pillows for pH electrode calibration. Each powder pillow prepares 50 mL of fresh buffer solution. Colour-coded for easy identification (Yellow). pH is 7.00 $\pm$ 0.02 at 25°C. Pack of 250 Powder Pillows

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Barva kodirana?:	Da
Koncentracija:	7,00 pH $\pm$ 0,02 pH pri 25 °C
Opis:	pH Buffer Calibration reagent
Package settings:	Powder Pillows
Parameter:	pH
pH vrednost(i):	7.00 $\pm$ 0.02 pH at 25°C
Pogoji skladiščenja:	10 °C - 30 °C
Posebni zapiski:	Discard prepared solution after use
Prostornina / velikost paketa:	250 pcs
Sledljivo po NIST-u:	Da