



Buffer Powder Pillows, pH 10.01 Colour-coded Blue, pk/50

Številka izdelka:

2227166

EUR Cena:

Kontakt

Na voljo

Nearno



pH Buffer standard solutions for laboratory and field analyses

pH Buffer Powder Pillows for pH electrode calibration. Each powder pillow prepares 50 mL of fresh buffer solution. Colour-coded for easy identification (Blue). pH is 10.01 $\pm$ 0.02 at 25°C. Pack of 50 Powder Pillows

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Barva kodirana?:	Da
Koncentracija:	10,01 $\pm$ 0,02 pH
Opis:	pH Buffer Calibration reagent
Package settings:	Powder Pillows
pH vrednost(i):	10.01 $\pm$ 0.02 pH at 25°C
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Discard prepared solution after use
Prostornina / velikost paketa:	50 pcs
Sledljivo po NIST-u:	Da