



Be Right™



Pufrska raztopina za pH; pH 10,01 (500 mL)

Številka izdelka: 2283649

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.

Predlagani nadomestki:

- 2283649.99

pH Buffer standard solutions for laboratory and field analyses

Kakovostna pufrska raztopina; označena z modro; sledljivost po NIST.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Barva kodirana?:	Da
Koncentracija:	10,01 pH ± 0,02 pH pri 25 °C
Opis:	Pufrska raztopina za pH; barvno označena za pH 10,01 (modra)
Package settings:	Bottle
Parameter:	pH
pH vrednost(i):	10.01 ±0.02 pH at 25°C
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Prostornina / velikost paketa:	500 mL
Sledljivo po NIST-u:	Da