



Be Right™



Standardna raztopina natrijevega klorida, 180 μS/cm (NIST)

Številka izdelka: 2307542

EUR Cena: Kontakt

Na voljo

Conductivity standard solutions for laboratory and field analyses

85,47 mg/L NaCl. Za elektrode za prevodnost

Primeren tudi za standardizacijo merilnikov prevodnosti; 180 μS/cm (μmho/cm). 100 mL

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Formula:	NaCl
Koncentracija:	NA
Opis:	Conductivity calibration solution
Package settings:	Bottle
Parameter:	Prevodnost
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	100 mL 100 mL
Shelf life:	24 months
Sledljivo po NIST-u:	Da
Volume / package size e Post Text:	100 mL